

REVISION HISTORY

Ī	Version	Author	Date	Revision Notes
	1.0	Nand Kishore		Initial Version

RS No: NA	System Description: PV Intake	Environment: (check one) Dev Test/Val Prod Prod
Test Script No:	Test Case / Name: PV Intake - AI Engine	Test Objective: Install and configure PV-Intake-Al Engine

Requirements Reference:	NA NA
Acceptance Criteria:	The objective for test script successfully met.

Author Approvals & Signatures

Title/Company Name

Name

Signature

Date

Author
(RxLogix)



Pre-Approvals			
Title/Company Name	Name	Signature	Date
SME/ Val Lead Reviewer (RxLogix)			
Quality (optional)* (RxLogix)	NA	NA	NA



SETUP DATA:

	Preparatory Work	Value	Comments / Notes	Verified By: Initial/Date
	<server ip=""></server>		Server IP on which user is performing installation	
	<port></port>	9999	Port on which user is performing installation	
	<ppk file=""></ppk>			
Test Data Input Table: Parameter Values	<pass key="" phrase=""></pass>		Pass phrase Key being used by user to log in to server using ssh command.	
values	<userid></userid>			
	<password></password>			
	<ami id=""></ami>		Image id to be used to create AWS servers	ami-074e2d6769f445be5



PREREQUISITES:

- 1. Tester is compliant and trainings have been completed which are required for the script execution.
- 2. Tester is able to login to the server using any ssh client and should have valid user id with sudo permissions on the server.
- 3. Must be connected to LAN network to avoid network issue.
- 4. Inbound and outbound traffic should be opened on <PORT> from aws security group, on which installation is being done
- 5. Tester is well-acquainted with the system.
- 6. Application server must meet the below hardware and operating system requirements, please note the below mentioned requirement may vary from client to client based on user load, usage and requirement.
- 7. Latest **docker_ai.zip** folders should be present in home directory.
- 8. Tester should take all file from S3 bucket.

System	Туре	Operating System	vCPU Cores	RAM	HDD	AWS Instance Type
PV Intake - AI Engine Production	Deployment Per Client Basis		4	16 GB or more	80GB	

PVI-4.8-IQ-004-1.0 Confidential Page 4 of 9



PROCEDURE

No.	Procedure	Expected Results	Actual Results	Pass / Fail	Verified By: Initial/Date
1.	Refer Test Data Input Table: Parameter Values for parameter values required in this step.	No error found, command executed	Attachment -		
	If you are using a Linux based system:	successfully.	Page		
	Log into the PV-Intake-Al Engine by executing the following command using the private key as mentioned in the prerequisites:				
	ssh -i <private key=""> <urerid>@ <server ip=""></server></urerid></private>				
	Note – Private key is the access key that will be used for logging into client's server.				
	Otherwise, if you are using a Windows system then log into <server ip=""> using a SSH client like Putty.</server>				
	Take multiple screenshots for all the commands, if required				



			Fail	Initial/Date
values required in this step.	No error found, command executed successfully.	Attachment - Page -		



No.	Procedure	Expected Results	Actual Results	Pass / Fail	Verified By: Initial/Date
3.	Refer Test Data Input Table: Parameter Values for parameter values required in this step.	Status should be online.	Attachment -		
	Step to check if API is up and running inside docker by executing below mentioned command:		Page		
	 unzip docker_ai.zip sudo docker build -t ai-engine:1.0 /home/<user_name>/docker_ai</user_name> sudo docker run -d -p 9999:9999restart alwaysname "ai" ai-engine:1.0 				
	Take multiple screenshots for all the commands.				



COMMENTS/REVIEW

TESTER COMMENTS		
Were all results Acceptable and test objective met?		
Yes No (If discrepancies were observed, refer to the Discrepancy Report Form(s) identified below)	Test Evidence Supporting documentation attached (# of attachments or N/A if not applicable).	Number of Attachments
Tester Comments:		
Completed By:		
(Signature/Date)		



Post Approvals			
Title/Company Name	Name	Signature	Date
SME/ Val Lead Reviewer (RxLogix)			
Quality (optional)* (RxLogix)	NA	NA	NA